LigaSure Advance™

Monopolar Tip Laparoscopic Sealer/Divider,
Pistol Grip

- en** Instructions for Use
- fr** Mode d'emploi
- pt** Instruções de uso
- de** Gebrauchsanleitung
- es** Instrucciones de uso
- it** Istruzioni per l'uso
- nl** Gebruiksaanwijzing
- sv** Bruksanvisning
- ru** Инструкция по применению
- zh** 使用说明

LigaSure Advance™

REF LF5544 **Monopolar Tip
Laparoscopic
Sealer/Divider,
Pistol Grip
5 mm shaft diameter,
44 cm shaft length**

For use with the Covidien ForceTriad™ energy platform with software version 3.2 or higher

For use with a maximum voltage of 288 Volts (peak) in the sealing (tissue fusion) mode

For use with a maximum voltage of 1,200 Volts (peak) in the dissection mode

Do not use if package is opened or damaged.

Defibrillation-proof type CF applied part

General Warnings

Warning

The product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection or product failure risks to the patient.

This device is designed for single use only. Subjecting it to reprocessing and/or multiple procedures can result in device failure. Failures can lead to compromised clinical outcomes such as seal failures and/or bleeding. Failure modes can include but are not limited to:

- Breakage during use
- Knife blade dullness or damage
- Sealing-jaw damage
- Electrical insulation failure resulting in electrical shocks, burns, sparks, and/or fire
- Intermittent or total device failure during use

This device is not designed for reprocessing and resterilization. Device features such as crevices, hinges, and lumens do not allow proper cleaning for sterilization, resulting in an increased risk of cross-infection, or exposure to disinfectant residues, endotoxins, or prions.

Warning

The LigaSure Advance LF5544 Monopolar Tip Laparoscopic Sealer/Divider is intended for use ONLY with the Covidien ForceTriad energy platform. Use of this instrument with other Covidien generators or with generators produced by other manufacturers may not result in the desired tissue effect, may result in injury to the patient or surgical team, or may cause damage to the instrument.

Do not use the LigaSure™ system unless properly trained to use it in the specific procedure being undertaken. Use of this equipment without training may result in serious unintended patient injury.

LigaSure vessel sealing systems have not been shown to be effective for tubal sterilization or tubal coagulation for sterilization procedures. Do not use these systems for these procedures.

Use caution during surgical procedures in which patients exhibit certain types of vascular pathology (atherosclerosis, aneurysmal vessels, etc.). For best results, apply the seal to unaffected vasculature.

Do not use this instrument on vessels in excess of 7 mm in diameter.

Confirm proper ForceTriad energy platform settings before proceeding with surgery.

Contact between an active instrument electrode and any metal object (hemostats, staples, clips, retractors, etc.) may increase the current flow and can result in unintended surgical effects such as an effect at an unintended site or insufficient energy disposition.

Electric Shock Hazard Do not connect wet accessories to the ForceTriad energy platform.

Fire Hazard Do not place instruments near or in contact with flammable materials (such as gauze or surgical drapes). Instruments that are activated or hot from use may cause a fire. When not using instruments, place them in a clean, dry, highly visible area, and not in contact with the patient. Inadvertent contact with the patient may result in burns.

The external surfaces of the instrument jaws may remain hot enough to cause burns after the radio frequency (RF) current is deactivated.

Do not use Valleylab™ mode for sealing. Doing so will result in an inadequate seal area.

Warning

Do not insert the instrument directly into a hand-assist port. Damage to the flexible insulation may occur.

Conductive fluids (e.g., blood or saline) in direct contact with or in close proximity to the instrument may carry electrical current or heat, which may cause unintended burns to the patient. Remove fluid from around the instrument jaws before activating the instrument.

For laparoscopic procedures, be alert to these potential hazards:

- Inadvertent activation or movement of the activated instrument outside of the field of vision may result in injury to the patient.
- Do not activate the instrument while instrument jaws are in contact with, or in close proximity to other instruments, including metal cannulas, as localized burns to the patient or physician may occur.
- Do not use hybrid cannulas that are comprised of both metal and plastic components. Capacitive coupling of RF current may cause unintended burns.
- Do not activate the instrument in an open circuit condition. Activate the system only when the instrument jaws are in full view and in direct contact with the target tissue to lessen the possibility of unintended burns.
- Carefully insert and withdraw instruments from cannulas to avoid possible damage to the devices and/or injury to the patient.
- In order to prevent damage to the flexible insulation proximal to the jaws, confirm the handles are in the fully closed position prior to insertion and extraction from the cannula.
- Use the appropriately sized cannula to allow for easy insertion and extraction of the instrument. Failure to do so may impact the integrity of the flexible insulation **(a)**.
- Cannulas with hard, non-beveled openings may cause the flexible insulation **(a)** to retract, which may compromise the insulation. If retraction occurs, the instrument must be discarded.

Always keep the external surface of the instrument jaws away from adjacent tissue while activating the instrument.

Warning

During a seal cycle, energy is applied to the area between the jaws. This energy may cause water to be converted to steam. The thermal energy of the steam may cause unintended injury to tissue in close proximity to the jaws. Care should be taken in surgical procedures occurring in confined spaces due to the possibility of steam formation.

Inspect the instrument and cords for breaks, cracks, nicks, or other damage before use. Failure to observe this caution may result in injury or electrical shock to the patient or surgical team, or cause damage to the instrument. If damaged, do not use.

Position instrument cords to avoid contact with the patient or other cords or cables. Do not wrap cords around metal objects. This may induce currents that could lead to shocks, fires, or injury to the patient or surgical team.

Keep the cords free from the jaw, shaft, and body of the instrument.

Examine all connections to the energy platform and all instruments and accessories before using. Improper connection may result in arcs, sparks, accessory malfunction, or unintended surgical effects.

Caution


Pediatric applications and/or procedures performed on small anatomic structures may require reduced power settings. The higher the current flow and the longer the current is applied, the greater the possibility of unintended thermal damage to tissue, especially during use on small appendages.

Before use, examine the energy platform and accessories for defects. Do not use cables or accessories with damaged insulation or connectors.

Intended Use


Each LigaSure Advance Monopolar Tip Laparoscopic Sealer/Divider instrument is for use in a single laparoscopic surgical procedure requiring monopolar dissection through tissue planes, creation of otomies (i.e., enterotomy or gastrotomy) and/or vessel sealing (tissue fusion) of vessels, pulmonary vasculature, tissue bundles and lymphatics. It is for exclusive use with the ForceTriad energy platform.

LigaSure Advance Instrument


- Ⓐ Flexible insulation
- Ⓑ Monopolar tip (dissecting electrode)
- Ⓒ Instrument jaws
- Ⓓ Rotation wheel
- Ⓔ LigaSure activation button (purple)
- Ⓕ Valleylab mode activation buttons (clear), one on the left side and one on the right
- Ⓖ Cutting trigger
- Ⓗ ElectroSurgery Smart™ connector (blue)
- Ⓘ LigaSmart™ connector (purple)
- ⓵ Handle


LigaSure Advance Connection Diagram on the ForceTriad Energy Platform

- Ⓗ ElectroSurgery Smart connector (blue)
- Ⓘ LigaSmart connector (purple)
- Ⓚ Monopolar 2 instrument receptacle ONLY
- Ⓛ Monopolar 2 screen
- Ⓜ LigaSure instrument receptacles
- Ⓝ LigaSure screen
- Ⓞ REM patient return electrode receptacle
- Ⓟ REM patient return electrode

①

②

- ① Do not use with the LigaSure vessel sealing system
- ② For use with the ForceTriad energy platform only

Product Description

The LigaSure Advance Monopolar Tip Laparoscopic Sealer/Divider is a specialty instrument that combines an advanced laparoscopic vessel sealing and division function with a monopolar tip in one instrument. The vessel sealing portion of the instrument permanently seals vessels up to and including 7 mm in diameter and tissue bundles without dissection or isolation.

Two output modes are selected on the handset with the following buttons:

- The purple LigaSure activation button **(e)** enables a bipolar LigaSure sealing function. See *LigaSure Function During Surgery* on page 5.
- The clear Valleylab mode activation buttons **(f)** enable a monopolar hemostasis with dissection function. See *Monopolar Function During Surgery* on page 6.

The instrument can also be used for tissue manipulation and blunt dissection. See *Tissue Manipulation and Blunt Dissection* on page 5.

Refer to the *ForceTriad Energy Platform User's Guide* for details of energy platform operation.

Setup

1. With the bar code on the LigaSmart connector **(i)** facing up, firmly insert it into either one of the LigaSure receptacles **(m)** under the LigaSure touchscreen **(n)** on the energy platform front panel.
2. With the bar code on the electro-surgery Smart connector **(h)** facing up, firmly insert it into the Monopolar 2 receptacle **(k)** under the center electro-surgery touchscreen on the energy platform front panel.

The ForceTriad energy platform recognizes the instrument type by reading the instrument's Smart connectors and sets the power output accordingly:

- The energy platform detects the LigaSmart connector **(i)** it sets the intensity setting to 2 bars in the display. If you have entered settings in the LigaSure touchscreen prior to connecting the LigaSure instrument, these settings will be reset to 2 bars.
- The energy platform also detects the electro-surgery Smart connector **(h)** and sets the power output to 15 watts in the Valleylab mode.

Notice

Monopolar function requires proper placement of a REM patient return electrode **(p)**. Refer to the patient return electrode instructions for use.

Tissue Manipulation and Blunt Dissection

The instrument can be used for manual dissection or to manipulate tissue with the jaws either open or closed. Rotating the instrument jaws provides another way to manipulate tissue.

Rotating the instrument jaws

With the instrument handle in the unlatched (open) position, the jaws can rotate 335 degrees.

Notice

Do not turn the rotation wheel when the handle is fully closed and latched. Product damage may occur.

1. Open the jaws by extending the handle **(j)** forward.

Note: If the handle is latched, squeeze it until it unlatches so you can then push the handle forward.

2. Turn the gray rotation wheel **(d)** on the instrument until the jaws are in the required position.

LigaSure Function During Surgery

LigaSure function can be used for these procedures:

- Sealing vessels and tissue bundles
- Cutting tissue

Warning

Place the vessel or tissue in the center of the jaws. To avoid incomplete sealing, do not place tissue in the jaw hinge or grasp tissue beyond the black line that indicates the rear edge of the electrode surface.

Do not activate the LigaSure function until the handle has been properly latched. Activating the function before this is done may result in improper sealing and may increase thermal spread to tissue outside the surgical area.

To assure proper function while sealing and cutting, eliminate tension on the tissue just prior to sealing and cutting.

Use LigaSure mode for sealing. Using Valleylab mode will result in an inadequate seal area and may damage the instrument.

Sealing Vessels And Tissue Bundles

1. Open the jaws by extending the handle **(j)** forward.

Note: If the handle is latched, squeeze it until it unlatches then push the handle forward.

- Grasp the intended vessel and/or tissue in the center of the jaws between the black lines indicating optimal tissue placement.

Warning

Do not attempt to seal over clips, staples, or other suture material as incomplete seals may be formed.

- Fully close the handle until it latches.
- Activate the instrument by either of the following methods:
 - Press and hold the purple LigaSure activation button (**e**) on the back of the instrument.
 - Depress and hold the corresponding round purple or orange pedal on the footswitch.

Note: Ensure that the footswitch you are using is connected to the same color port as the instrument (purple to purple, orange to orange).

A continuous tone sounds to indicate that the vessel or tissue is being sealed. When the activation cycle is complete, a two-pulsed end tone sounds and the energy platform discontinues RF output.

- Release the purple LigaSure activation button on the instrument or purple or orange foot pedal when the sealing cycle is complete.
- To seal adjacent tissue, overlap the edge of the existing seal.

See *Cleaning the Instrument During Use* on page 7 for cleaning instructions.

Cutting Tissue

Caution

Energy based devices, such as electro-surgical pencils or ultrasonic scalpels, that are associated with thermal spread should not be used to transect seals.

Do not overfill the jaws of the instrument with tissue. This may compromise the cutting function.

Confirm that the jaws have reached the closed position and are locked (the handle is latched) before activating the cutter.

Notice

If the instrument shaft is bent, the knife will not engage. Discard and replace the bent instrument.

Do not engage the cutting mechanism over clips, staples, suture material, or other foreign objects as damage to the cutter may occur.

Do not apply force to the shaft of the instrument causing tension or bowing as this could make the knife difficult to deploy and the trigger may not return to its normal position.

- Prior to cutting the seal, the surgeon may inspect the vessel or tissue to ensure proper sealing.
- Verify that the handle (**j**) is fully closed and latched.
- To engage the cutting mechanism, gently pull the cutting trigger (**g**) completely back toward the rear of the instrument.
- Release the cutting trigger to return to the starting position.
- Open the jaws by squeezing the handle until it unlatches, then push the handle completely forward from the body of the instrument.

See *Cleaning the Instrument During Use* on page 7 for cleaning instructions.

Monopolar Function During Surgery

Monopolar function enables hemostasis with dissection and requires a REM™ patient return electrode (**p**).

Warning

Tissue buildup (eschar) on the tip of an active electrode may create embers that pose a fire hazard, especially in oxygen enriched environments. Keep the electrode clean and free of all debris. This will reduce the need for additional energy and decrease thermal spread.

Prior to activating electro-surgical current, verify that the flexible insulation completely covers the metal portion of the shaft directly proximal to the jaws to reduce the potential for unintended tissue effect.

Activate the instrument **only** when ready to deliver electro-surgical current and when the monopolar tip is in view (especially if looking through an endoscope).

Warning

Do not activate the monopolar tip while grasping tissue in the jaws. Doing so will result in unintended tissue effect.

Confirm proper electrosurgical settings prior to and during a procedure. Use the lowest power settings that achieve the desired effect. If increased power settings are requested, check the patient return electrode and all accessory connections before major power setting adjustments.

Conductive fluids (e.g., blood or saline) in direct contact with an active electrode or in close proximity to an active accessory may carry electrical current and cause unintended burns to the patient. This can happen as a result of either direct coupling with the active electrode or capacitive coupling between the active electrode and the external surface of the electrode insulation. To prevent unintended burns in the presence of conductive fluids:

- Always keep the external surface of the active electrode away from adjacent tissue while activating the instrument.
- Clear conductive fluid from the electrode before activating the electrosurgical device.

Activating the monopolar tip for extended periods of time may damage the instrument jaws and cause the jaws to retain heat.

Always use the lowest power setting that achieves the desired surgical effect. Use the active electrode for the minimum time necessary in order to reduce the possibility of unintended burn injury.

Deactivate the electrosurgical unit **before** the tip leaves the surgical site.

Danger: Explosion Hazard Do not use electrosurgery in the presence of flammable anesthetics.

Fire Hazard Sparking and heating associated with electrosurgery can be an ignition source. Keep gauze and sponges wet. Keep electrosurgical electrodes away from flammable materials and oxygen (O₂) enriched environments.

Use of electrosurgery in O₂ enriched environments increases the risk of fire. Therefore, take measures to reduce the O₂ concentration at the surgical site.

Avoid enriched O₂ and nitrous oxide (N₂O) atmospheres near the surgical site. Both O₂ and N₂O support combustion and may result in fires and burns to patients or surgical personnel.

Warning

Prevent pooling of flammable fluids and the accumulation of flammable or oxidizing gases or vapors under surgical drapes or near the surgical site.

Do not activate the instrument until flammable vapors from skin prep solutions and tinctures have dissipated.

Avoid the accumulation of naturally occurring flammable gases that may accumulate in body cavities such as the bowel.

Activating Monopolar Function

The power in Valleylab mode is controlled via the Monopolar 2 screen **(I)** on the generator. Changes can be made in 1 watt increments. See the *ForceTriad Energy Platform User's Guide* for information about using the touchscreen.

Note: The instrument default setting delivers 15 watts in Valleylab mode.

Notice

Before activating the monopolar electrosurgical function, a REM patient return electrode **(p)** must be properly placed on the patient. If a REM alert occurs on the ForceTriad energy platform, refer to the patient return electrode instructions for use for proper placement and usage.

1. Verify that the instrument handle **(j)** is in the fully closed and latched position.
2. Activate the monopolar tip by pressing one of the clear Valleylab-mode activation buttons **(f)** on either side of the LigaSure Advance instrument. A tone will sound indicating that energy is being delivered.

See *Cleaning the Instrument During Use* on page 7 for cleaning instructions.

Important

Monopolar electrosurgery cannot be activated using a footswitch.

Cleaning the Instrument During Use**Warning**

Do not activate the instrument or the cutting trigger while cleaning the jaws or monopolar tip. Injury to operating personnel may result.

Important

Keep the instrument jaws and monopolar tip clean. Buildup of eschar may reduce sealing, cutting, and/or dissection effectiveness.

- Gently wipe the jaws, monopolar tip, and inside surfaces with a wet gauze pad.
- Do not attempt to clean the flexible insulation. Cleaning may damage insulation.
- Do not clean the instrument jaws or monopolar tip with a scratch pad or other abrasive material.

Troubleshooting

Methods for troubleshooting some situations appear here. For situations not covered here, see the *ForceTriad Energy Platform User's Guide*.

Alert Situations

Note: If both the handswitch and footswitch pedal are activated during the same seal cycle, the ForceTriad energy platform delivers energy from the activation source it detects first.

When an alert condition occurs, energy delivery is terminated, but will be available immediately after the alert condition has been corrected.

Check Instrument

A *Check Instrument* alert condition produces six pulsed tones requiring the user to check the instrument and the LigaSure touchscreen on the ForceTriad energy platform displays a message.

If the *Check Instrument* message appears on the LigaSure touchscreen, the user should:

1. Release the footswitch pedal or LigaSure activation button.
2. Open the instrument jaws and inspect for a successful seal.
3. Follow the suggested remedies on the *Check Instrument* screen.

If possible, reposition the instrument and regrasp tissue in another location, then reactivate the seal cycle.

Reasons for Check Instrument Alert

Regrasp thicker tissue – Indicates the user is grasping thin tissue or not enough tissue; open the jaws and confirm that a sufficient amount of tissue is inside the jaws. If necessary, increase the amount of tissue and repeat the procedure.

Check for clips / Regrasp tissue – Avoid grasping objects, such as staples, clips or encapsulated sutures in the jaws of the instrument.

Clean jaws and tips – Use a wet gauze pad to clean surfaces and edges of instrument jaws.

Remove excess fluids – Pooled fluids around the instrument tip; minimize or remove excess fluids.

Endpoint Not Reached

The *Endpoint Not Reached* alert condition produces four pulsed tones. Also, the LigaSure touchscreen on the ForceTriad energy platform displays an alert message that instructs the user about corrective actions.

If the *Endpoint Not Reached* message appears on the LigaSure touchscreen, the user should:

1. Release the footswitch pedal or activation button.
2. Reactivate the seal cycle without removing or repositioning the jaws from their original position.

Reasons for Endpoint Not Reached Alert

Maximum seal cycle time has been reached – The system needs more time and energy to complete the seal.

Seal cycle was interrupted before the cycle was complete – The footswitch or LigaSure activation button was released before the seal cycle was complete.

After Surgery

1. Turn off the ForceTriad energy platform.
2. Disconnect the Electrosurgery Smart connector (**h**) and the LigaSmart connector (**i**) from the ForceTriad energy platform.
3. Discard the LigaSure Advance sealer/divider instrument using standard precautions and according to the facility procedures. The instrument is not able to withstand resterilization. **Do not resterilize.**

Warning

Do not reuse or resterilize the LigaSure Advance sealer/divider instrument or other accessories labeled “disposable” or “single use only.” See the detailed warning on page 2.

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© 2010 Covidien.
COVIDIEN, COVIDIEN with logo and the Covidien logo are U.S. and/or internationally registered trademarks of Covidien AG.
Other brands are trademarks of a Covidien company.
Tyco Healthcare UK Ltd., Gosport, PO13 0AS, UK.